



CL 701



CL 702

ANTI BIOTIC ZONE READER CL 701 & CL 702

Specifications:

- This is very compact and designed for measuring the diameter of finished zone in Petri dish, for determining the strength of antibiotic materials.
- This Antibiotic Zone Reader measures the diameter of the inhibited zone within the range of 0 to 35 mm with accuracy of 0.2 mm.
- The light from inside source if the base passes up through the transparent and semitransparent of the Agar then to a reflecting mirror situated in an arm above the unit.
- The mirror reflects the light to a glass prism mounted at the magnified image of the zone of inhibitions is clearly invisible on the prism.

CROMTECH INDIA

1st Floor Begum Niwas, Opp Mercantile, Bank Society Colony, 5th Natwar Nagar Road, Jogeshwari (East), Mumbai 400 060.

Maharashtra - India.

Mobile: +91 98693 82134 / 81088 13320

Telefax: +91 22 2838 6028

Email: cromamumbai@mail.com/ cromtechindia@mail.com/ info@cromtechindia.com/ sales@cromtechindia.com

Website: www.cromtechindia.com